


Catalog No. M-1100




Patent Applied

# R & F<sup>®</sup> *Listeria* sp./*Listeria monocytogenes* Chromogenic Plating Medium


A Selective & Differential Medium For Detecting  
Nonpathogenic & Pathogenic *Listeria* On a Single Plate

Presumptively positive colonies of nonpathogenic *Listeria* sp. appear pink, whereas *Listeria monocytogenes* and *Listeria ivanovii* appear as blue-green to blue-violet colonies with or without a precipitate after 42-48 hours at 35°C.



*Listeria*  
*monocytogenes*

*Listeria innocua*



Fluorescence

Acid production  
from rhamnose

R & F<sup>®</sup> *Listeria* sp./*Listeria monocytogenes*  
Chromogenic Plating Medium Powder

R & F<sup>®</sup> *Listeria monocytogenes*  
Confirmatory Medium  
Separates *L. monocytogenes* from *L. ivanovii*

**R & F® *Listeria* sp./*Listeria monocytogenes*  
Chromogenic Plating Medium at 35°C for 42-48 hours**

Bacterial Species	# of Strains	Colonial Morphology
<i>Listeria monocytogenes</i>	39	Convex, 1-2 mm, blue-green to blue-violet, ± precipitate
<i>Listeria ivanovii</i>	4	Convex, 1-1.5 mm, dark blue-green, large ppt.
<i>Listeria innocua</i>	6	Convex, 1-2 mm, pink, no ppt.
<i>Listeria welshimeri</i>	2	Convex, 1-2 mm, pink, no ppt.
<i>Listeria seeligeri</i>	1	Convex, 1-2 mm, pink, no ppt.
<i>Listeria grayi</i>	1	Convex, 1-2 mm, pink, no ppt.
<i>Bacillus cereus</i> / <i>B. thuringiensis</i>	3	No Growth
<i>Enterococcus</i> sp.*		Scant Growth; pinpoint; clear
Gram positive sp.**		No Growth
Gram negative sp.***		No Growth
Yeasts sp.****	3	No Growth

\* *Enterococcus faecalis*, *E. faecium*, and *E. avium*

\*\* Includes: *Bacillus circulans*, and *B. subtilis*; *Staphylococcus aureus*, *S. epidermidis*, and *S. saprophyticus*; *Lactobacillus acidophilus* and *L. plantarum*; and *Pediococcus cerevisiae*

\*\*\* Includes: *E. coli* (2 strains) and *E. coli* 0157:H7 (1 strain); *Enterobacter aerogenes*; *Citrobacter freundii*; *Shigella sonnei*; *Morganella morganii*; *Providencia alcalifaciens*; *Pantoea agglomerans*;

*Enterobacter sakazakii*, *Klebsiella pneumoniae* and *K. ozaenae*

\*\*\*\* *Zygosaccharomyces bailii* and *Z. rouxii*; *Candida albicans*

**Advantages of R & F® *Listeria* sp./*Listeria monocytogenes* Chromogenic Detection System**

- This highly differential system is based on specific chromogenic substrates in the plating medium that simultaneously differentiate both nonpathogenic *Listeria* species (*L. innocua*, *L. seeligeri*, *L. welshimeri*, and *L. grayi*) as pink colonies and the pathogenic *Listeria* species (*L. monocytogenes* and *L. ivanovii*) as blue-green to blue-violet colonies on a single plate in 42-48 h at 35°C.
- Unlike selective/differential agar for *Listeria* species that depend on the detection of only β-glucosidase activity and produce a single color for colonies for all *Listeria* species, the differentiation mechanism in this system works by the combination of indoxyl derivative chromogenic substrates that produce positive color reactions for colonies of nonpathogenic *Listeria* sp. that are pink due to their β-glucosidase activity, and blue-green to blue-violet for the pathogenic species depending on the strain-specific balance of β-glucosidase (pink) and phosphatidylinositol-specific-phospholipase C (blue) activities.
- The agar surface of the plates have an opaque white background that facilitates differentiation of colored colonies of both groups of *Listeria* organisms growing on the surface of the plate.
- A rapid and convenient fluorogenic test and acid from rhamnose are provided in the system that differentiates *Listeria monocytogenes* from *Listeria ivanovii* within 6 hours.

**REFERENCES**

1. Restaino, L., E.W. Frampton, W.C. Lionberg, and A.L. Restaino. 2006. A multi-chromogenic agar for the dual detection of nonpathogenic and pathogenic *Listeria* species. Poster No. P3-64, IAFP 2006.
2. Restaino, L., E.W. Frampton, W.C. Lionberg, and A.L. Restaino. 2007. Detection of non pathogenic and pathogenic *Listeria* species by use of a chromogenic agar. *Food Prot. Trends*. 27:592-596.
3. Swiech, R. et al. 2008. Efficacy of a chromogenic plating medium for detaching *Listeria* species from environmental samples. Poster No. P3-16, IAFP 2008.

**ORDERING INFORMATION:**

**M-1100 R & F® *Listeria* sp./*Listeria monocytogenes* Chromogenic Plating Medium**

**M-1110 R & F® *Listeria* sp./*Listeria monocytogenes* Supplement for Plating Medium**

**M-1150 R & F® *Listeria* sp./*Listeria monocytogenes* Chromogenic Detection System**

**M-0520 R & F® *Listeria monocytogenes* Confirmatory Medium**

R & F Products, Inc.  
2725 Curtiss Street  
Downers Grove, IL 60515

Phone: 630.969.5300  
Fax: 630.969.5303  
E-mail: rf\_labs@sbcglobal.net  
Web Site: www.rf-labs.com